



# HYDRO-GAZ-MED


PRODUCENT WYROBÓW MEDYCZNYCH

ul. Willowa 40, 05-205 Dobczyn

☎/☎ 22 787 65 60 / ☎ 606 348 532

🌐 www.gazmed.pl / ✉ gazmed@gazmed.pl

## Gas analyzer *information leaflet*



Purpose of this device is to continuously monitor parameters of medical compressed air and medical oxygen according to EN ISO 7396-1:2016, European Pharmacopeia and HTM 02-01.

### Technical data:

- required gas flow: 5-10 l/min
- inlet pressure: 5bar
- dimensions: 320x400x150 (mm) (WxHxD)
- weight: ~5kg
- power supply: 230VAC

### Measuring ranges:

- Oxygen purity 0-100% V/V
- Pressure 0-16bar
- Dew point -70-100°C
- Flow 0-800 l/min
- H<sub>2</sub>O 20-1000 ppm
- CO<sub>2</sub> 0-2000 ppm
- NO<sub>2</sub> 0-500 ppm
- CO 0-2000 ppm
- NO 0-2500 ppm
- SO<sub>2</sub> 0-500 ppm