



Switchover system for cylinder supply "PNEUMAT I"

PNEUMAT I as 2 cylinder banks supply system



Dimensions
400x700x150
(WxHxD)



Remote alarm unit

- Capacity: max. 30 m³/h
Inlet pressure: 200 bar (2900 psi)
Outlet pressure: 4-6 bar (58-87 psi)
Automatic change-over (pneumatically)
Connection to BMS and MODBUS RTU
Working modes:
- ▶ 2 cylinder banks
 - ▶ 2 cylinder banks + 1 reserve bank
 - ▶ Concentrator + 2 cylinder banks
 - ▶ Tank with liquid gas + 2 cylinder banks
 - ▶ Air compressor + 2 cylinder banks
- Class IIb medical device.