



Switchover system for cylinder supply "PNEUMAT II"

PNEUMAT II as 2 cylinder banks + reserve bank supply system



Dimensions
450x800x200
(WxHxD)



Remote alarm unit

- Capacity: max. 50 m³/h
 Inlet pressure: 200 bar (2900 psi)
 Outlet pressure: 4-6 bar (58-87 psi)
 Automatic change-over (electronically)
 Connection to BMS and MODBUS RTU
 Working modes:
- ▶ 2 cylinder banks
 - ▶ 2 cylinder banks + 1 reserve bank
 - ▶ Concentrator + 2 cylinder banks
 - ▶ Tank with liquid gas + 2 cylinder banks
 - ▶ Air compressor + 2 cylinder banks
- Class IIb medical device.