


High pressure manifold system

Technical data:

- ▶ up to 10 cylinders in single manifold
- ▶ shut off valve
- ▶ exhaust valve
- ▶ non-return valves for each cylinder
- ▶ gas specific connections




Gas coding colours:




White - Oxygen



Black - Compressed air




Grey - Carbon Dioxide



Blue - Nitrous Oxide

Manifolds delivered along with headers.



2 cylinder header



Single cylinder header

Technical data:

- ▶ maximum inlet pressure 200 bar (2900 psi)
- ▶ easily adjustable pigtail shape
- ▶ inlet/outlet connections according to gas coding standards

Technical data:

- ▶ maximum inlet pressure 200 bar (2900 psi)
- ▶ length ~80cm
- ▶ inlet/outlet connections according to gas coding standards