

Area valve service unit



Technical data:

- ▶ HTM 02-01 compliant
- ▶ pressure sensors for each medium
- ▶ shut-off valve with lock for each medium
- ▶ emergency NIST connectors either side of the valve
- ▶ gauge for each medium
- ▶ physical separation
- ▶ front panel open/close detection
- ▶ alarm unit with 7" LCD touch screen
- ▶ connection to BMS and MODBUS RTU

Working pressure:

- ▶ medical gases: 0-10 bar (0-145 psi)
- ▶ vacuum: 0- -0.9 bar (0- -13 psi)

