

Gas outlet panel "SPG"



Gas outlet panel "SPG-5"



Gas outlet panel "SPG-6"



Gas outlet panel "SPG-7"



Gas outlet panel "SPG-9"

Gas outlet panels
(under plaster version)



Gas outlet panel (on-plaster version)

Medical gas outlets mounted in panel for under plaster installation are used to acquire medical gases and vacuum. Panel can be equipped with any number of gas outlets either with or without anesthetic gas scavenging system or AIRMOTOR.

Gas outlets can be in AGA SS 875 24 30, DIN 13260-2, NF S 90-116 and BS 5682 connection standard.

In our products we use gas outlets produced by renowned company GREGGERSEN GASETECHNIK.